




mindray
healthcare within reach

CAL 6000

Línea de análisis celular:
hematología y frotis automatizados
a su alcance

- Sistema modular de mesada.
- Combinación de dos contadores hematológicos BC-6000 o BC-6200 y una unidad SC-120, dedicada a la elaboración de frotis sanguíneos.
- Lector de código de barras para racks que distribuye las muestras según las necesidades del usuario.
- Velocidad: 220 hemogramas/hora con 120 frotis/hora.
- Recuento de eritroblastos en todos los hemogramas.
- Repetición y test reflejos automáticos.
- Análisis de líquidos biológicos sin reactivos adicionales.
- Software labXpert: gestiona todas las muestras desde una única computadora.



REPRESENTANTE EN ARGENTINA

GEMATEC
equipamiento para medicina



Int. Avalos 3651 | (1605) | Munro
Buenos Aires, Rep. Argentina
Tel./Fax: (54 11) 4512 5666
ventas@gematec.com.ar
www.gematec.com.ar

mindray

Mindray is listed on the NYSE under the symbol "MR"
Mindray Building, Keji 12th Road South,
High-tech Industrial Park,
Nanshan, Shenzhen 518057, P.R. China
Tel: +86 755 21888998 Fax: +86 755 26582680
E-mail: intl-market@mindray.com
Website: www.mindray.com



MINDRAY is a trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Specifications subject to changes without prior notice.
©2010 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.
P/N: SPA-BC2800-210285x4-20110310